



# VCX 130PB

130 Watts

Processing Volume: 200µl – 150mL\*



The **VCX 130PB** is a 130 watt system with pulse button control. The probe can be handheld and ultrasonic energy manually activated by depressing the button on the side of the converter.

This is a cost effective alternative to the **VCX 130** version. The 1/8" probe (0.5 – 15mL volumes) is included and a variety of accessories are available for various applications and sample volumes.

## FEATURES:

- Variable power output
- Pulse button control
- Manual pulse mode
- Wattage and Energy display



(stand is optional)

Ordering information: Model no. **VCX 130PB**  
130 watt ultrasonic processor with 1/8" diameter probe and tool kit.

## SPECIFICATIONS

Power Output	130 watt maximum	Converter Cable Length	Permanently attached 6 ft. (1.8m)
Frequency	20kHz		
Dimensions	H x W x D: 4.5 x 9.75 x 12.5" (115 x 250 x 320mm)		
Ultrasonic Converter	Part #CV188 (with pulse button)	Tool Kit	Part #888-00026 open end wrench (2 included)
Diameter	1.25" (32mm)		
Length	5.75" (146mm)		
Standard Probe	Part #630-0422	Optional:	
Tip Diameter	1/8" (3mm)	Stand with Clamp	Part #830-00460
Length	5.4" (137mm)	Sound Enclosure	Part #830-00451
Material	Ti6Al4V titanium alloy		

\*The actual processing volume of a specific probe or accessory is application specific.