



The VC 50 is a small and compact 50 watt ultrasonic processor. The probe can be held by hand and ultrasonic energy is activated by turning the amplitude knob on the front panel.

The 1/8" probe (0.5 - 15mL volumes) is included and a variety of accessories are available for various applications and sample volumes.

## **FEATURES:**

- Compact design
- Variable power output
- Low cost



(stand is optional)

Ordering information: Model no. VC 50

50 watt ultrasonic processor with 1/8" diameter probe and tool kit.

## **SPECIFICATIONS**

Power Output Frequency Dimensions	50 watt maximum 20 kHz H x W x D: 8 x 8 x 5" (203 x 203 x 127mm)	Converter Cable Length	Permanently attached 5 ft. (1.5m)
Ultrasonic Converter Diameter Length	Part #CV18 1.25" (32mm) 5.75" (146mm)	Tool Kit	Part #888-00026 open end wrench (2 included)
Standard Probe Tip Diameter Length Material	#630-0422 1/8" (3mm) 5.4" (137mm) Ti6Al4V titanium alloy	Optional: Stand with Clamp Sound Enclosure	Part #830-00460 Part #830-00451

<sup>\*</sup>The actual processing volume of a specific probe or accessory is application specific.