



The VC 505 is a cost effective alternative to the VCX model. This 500 watt system offers the same processing power and choice of accessories as the VCX500 without temperature monitoring or energy set point programming.

A standard ½" diameter probe with replaceable tip is included and many optional accessories are available to meet the needs of any application.



(stand is optional)

Ordering information: Model no. VC 505

500 watt ultrasonic processor with 1/2" diameter probe and tool kit.

SPECIFICATIONS

Power Output Frequency Dimensions	500 watt maximum 20kHz H x W x D: 9.25 x 7.5 x 13.5" (235 x 190 x 340mm)	Converter Cable Length	Part #201-0300 6 ft. (1.8m)
Ultrasonic Converter Diameter Length	Part #CV334 2.5" (64mm) 7.25" (183mm)	Tool Kit	Part #381-0005 Includes: 2 spanner wrenches and 1 open-end wrench
Standard Probe Tip Diameter Length Material	Part #630-0220 0.5" (13mm) with replaceable tip 5.5" (139mm) Ti6Al4V titanium alloy	Optional: Stand with Clamp Sound Enclosure	Part #830-00459 Part #830-00427

*The actual processing volume of a specific probe or accessory is application specific.

FEATURES:

- Programmable
- Variable power output
- 10 hour timer
- Pulse mode
- Wattage and Energy display