



The VCX 130 is a 130 watt, programmable ultrasonic processor and our most popular small volume system.

This model offers the same programmability as larger **VCX** versions.

The $\frac{1}{4}$ " probe (10 – 50ml volumes) is included and a variety of accessories are available for various applications and sample volumes.

FEATURES:

- Programmable
- Pulse mode
- 10 hour timer
- Variable power output
- Wattage and Energy display



Ordering information: Model no. VCX 130

130 watt ultrasonic processor with 1/4" diameter probe and tool kit.

SPECIFICATIONS

Power Output Frequency Dimensions	130 watt maximum 20kHz H x W x D: 4.5 x 9.75 x 12.5" (115 x 250 x 320mm)	Converter Cable Length	Permanently attached 5 ft. (1.5m)
Ultrasonic Converter Diameter Length	Part #CV18 1.25" (32mm) 5.75" (146mm)	Tool Kit	Part #888-00026 open end wrench (2 included)
Standard Probe Tip Diameter Length Material	Part #630-0435 1/4" (6mm) 5.4" (137mm) Ti6Al4V titanium alloy	Optional: Stand with Clamp Sound Enclosure	Part #830-00460 Part #830-00451

^{*}The actual processing volume of a specific probe or accessory is application specific.